REF 3031.1.00

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1 | Delivery

- 1 x Pelvic Muscle Trainer
- 1 x Instruction manual

Read the instructions for use completely and carefully!

2 | Intended purpose

The Pelvic Muscle Trainer

- serves to strengthen the pelvic floor muscles

Pelvic floor training increases the strength and tone of the muscles and makes it much easier to control the urge to urinate. Thus, targeted pelvic floor training is the best remedy - not only as a preventive measure, but also to combat existing urinary incontinence. The pelvic floor trainer is placed a few centimetres below the pelvic floor between the thighs, then the thigh muscles are squeezed. The squeezing is intended to strengthen the pelvic floor muscles.

Applicable to:

- Disc discomfort
- Prevention of uterine subsidence
- reduces or prevents urinary incontinence

For private use.

3 | Contraindications / Not suitable

- After recent operations, especially in the pelvic floor area, training should only be carried out after consultation with the attending physician. How long after an operation will the corresponding muscles remain. The doctor decides whether to take care of things.
- Allergic reactions to the material


Not suitable for children under 3 years!

4 | Notes on safe use of the product

- Side effects that affect the entire organism are not yet known. Proper use is assumed. Product liability is excluded in the event of improper use.
- People with limited physical, sensory or mental abilities should only use the pelvic floor exercise under supervision.
- If your symptoms worsen or you notice unusual changes, consult your doctor immediately.
- For additional questions consult your doctor!
- Due to the material composition, it can be worn directly on the skin as long as there are no external or internal injuries. However, it should only be used with clothing during training.
- To avoid pinching the skin, only use the pelvic floor exercise on clothing
- Examine the appliance frequently. If there is any damage visible If the supply cord is damaged, it must be replaced by the manufacturer or a qualified person, to avoid a hazard.

5 | Application


1. Warm up before and wear comfortable clothes while using Pelvic Muscle Trainer.
2. Place Pelvic Muscle Trainer between thighs, few inches below pelvic floor.
DO NOT place against pelvic floor.
3. Center the Pelvic Muscle Trainer with Wing flaps touching front or back of thighs depending which is more comfortable
4. Maintain proper posture. Stand up straight with feet close together, knees slightly bent.
5. While keeping knees, hips and shoulders parallel, close Wing Enclosure Area by squeezing Wing Flaps with upper thighs, buttocks and pelvic floor muscles.
6. For beginners, squeeze und hold the Pelvic Muscle Trainer for 3 seconds and then rest 7 seconds. Try to perform 6 repetitions.

You may find it difficult to close the Wing Enclosure Area at first. This is OK because you are still activating the PC muscles even when you get the Wing Flaps halfway closed.

7. Gradually increase usage and squeeze hold time. Start by using the Pelvic Muscle Trainer three times a week and increase to twice daily. Give with enough time for the muscles to recuperate. Avoid Straining Themuscles.
8. Experienced users should be able to perform 10 repetitions. Holding 25 seconds resting 5 seconds
9. As long as you maintain posture, exercise regularly and keep the Pelvic Muscle Trainer positioned properly the strength will come automatically.

Always make sure you close the Flex Flaps as far as you can with each repetition. The training is efficient and has no side effects. Effects are individually variable depending on the regularity of use.

Other applications for training the inner thigh muscles or the chest muscles:

6 | Cleaning

- Store the device in a dry, clean place
- After each use, clean the device with a dry or damp cloth

7 | Technical data

Model: SKE2

Material: Polypropylen

Weight: 205 gr

Dimensions: 20 cm x 23 cm x 8 cm

This product may be safely disposed of in normal home trash

8 | Warranty

The guarantee lasts for 24 months from the date of purchase for manufacturers' error in material and quality. If it is used for another purpose or operated incorrectly, no liability for possible damages can be accepted.

Address for the guarantee and service:

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Hydas GmbH
c/o atrikom fulfillment GmbH • Carl-Zeiss-Straße 10