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- 1× Foot massager (main unit)
- 1× Charging cable
- 1× User manual

## 1. Product Overview

This foot massager combines air compression massage, vibration massage, roller massage and an optional heating function. It is designed to provide relaxation, improve blood circulation and relieve foot fatigue.

The device can be used as a **foot massager** or as a **stool** when the top cover is closed.

### Main Components:

- Upper cover
- Lower body
- Foot chambers
- Control panel with LED indicators
- Power connector

## 2. Technical Data

### Massage Functions:

- Airbag massage
- Vibration massage
- Roller massage
- Heating function

### Programs:

- Green indicator: Airbags + Vibration + Rollers
- Yellow indicator: Airbags + Rollers
- Red indicator: Vibration

### Message Intensity (3 Levels):

- Green: Low intensity
- Yellow: Medium intensity
- Red: High intensity

### Heating Function:

- Switchable via the Heat button
- Provides gentle warmth during the massage

### Power Supply / Battery:

- Powered via adapter
- Two versions:
  1. Without battery → operation only with adapter
  2. With battery → cordless use after charging
- Charging time: approx. 3 hours

### Battery level indicator:

- Green: 70–100 %
- Yellow: 30–70 %
- Red: 0–30 %

### **3. Preparation Before Use**

#### **Installation Position**

- Keep at least 1 meter distance from TVs, radios or other electronic devices.

#### **Caution**

- Keep the adapter cable straight to avoid damage.
- Do not exceed the recommended charging time.
- Do not use non-original power adapters.

#### **Grounding**

- Device must be operated with proper grounding.
- Do not use on uneven surfaces.

#### **Environment**

- Keep away from moisture (e.g. bathroom).
- Avoid high temperatures.

### **4. Basic Operation**

#### **Power Supply**

- Version without battery → use with adapter only
- Version with battery → cordless operation possible

#### **Using as a Stool**

- Close the upper cover.
- The device can be used as a seat.

#### **Using as a Foot Massager**

- Open the cover and place your feet into the foot chambers.

### **5. Controls and Functions**

#### **Power Button**

- Turn device on/off
- LED battery indicator:
  - Green = 70–100 %
  - Yellow = 30–70 %
  - Red = 0–30 %

#### **Program Button**

- Green = Airbags + Vibration + Rollers
- Yellow = Airbags + Rollers
- Red = Vibration

### **Heat Button**

- Press once: Heating ON
- Press again: Heating OFF

### **Intensity Button**

- Switch between 3 intensity levels
  - Green = low
  - Yellow = medium
  - Red = high

## **6. Usage Notes**

- Recommended session duration: 15–20 minutes
- Allow at least 30 minutes between sessions
- Stop immediately if you feel pain or numbness
- Do not move your feet abruptly during use

## **7. Cleaning & Maintenance**

### **Cleaning**

- Disconnect device from power source
- Clean with a dry or slightly damp cloth
- Do not use aggressive cleaners

### **Maintenance**

If excessive dust accumulates, the unit may be cleaned with water or a neutral detergent.

Ensure all components are completely dry before reassembly.

### **Storage**

- Store in a dry and ventilated place
- Avoid heat and sunlight
- Do not bend or pinch the cable

## **8. Safety Instructions**

- Do not use in damp environments
- Do not use if the cable is damaged
- Do not open or repair the device
- Do not use while sleeping
- Not suitable for children under 14
- People with pacemakers, implants or health issues should consult a doctor

## **9. Disposal**

Electrical devices must not be disposed of in household waste.  
Dispose of according to local electronic waste regulations.

## **10. Warranty**

- Warranty according to legal regulations
- No warranty in case of:
  - Improper use
  - Unauthorized repairs
  - Water damage

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SHP Energy GmbH